

GOT PARTS?

TRIED AND TRUE SYRINGES.

*Press-Ring, A-dec style,
single button, Filter . . .*

We have them all.

(And we ship cookies with every order)

C-1900 "Forma" Press-ring Syringe - MSRP \$70

- ▶ Includes C-1623 indexing tip
- ▶ Trim, lightweight head
- ▶ Press-ring Quick-release tip
- ▶ Also accepts standard autoclavable tips
- ▶ Smooth lines enable thorough wipe down
- ▶ Cross pin allows easy button replacements
- ▶ Parts cross-compatible with domestic press-ring syringes



C-1915 "Forma" Single-button Syringe MSRP \$70

- ▶ Same great features as our C-1900 Air/Water Syringe with just one button!

C-1700 "Contour" Syringe MSRP \$70

- ▶ Trim, lightweight head
- ▶ Traditional push/pull tip release
- ▶ Accepts standard autoclavable tips (C-1624)
- ▶ Smooth lines enable thorough wipe down
- ▶ Cross pin allows easy button replacements



C-1715 "Contour" Filter Syringe MSRP \$90

- ▶ Single button offers variable air control
- ▶ Filters help remove air-borne moisture/contaminants
- ▶ Traditional push/pull tip release
- ▶ Filters commonly available (C-1716-10)
- ▶ Accepts all standard autoclavable tips (C-1624)
- ▶ Cross pin allows easy button replacements
- ▶ Includes 7 Feet of grey hookup tubing



Delivering cookies is great, but remember . . .

We deliver Parts



 **PARTS WAREHOUSE**
Let us be your dental resource

1-800-817-6704
www.dentalparts.com
CustomerService@Dentalparts.com
100% American Owned