

# GOT SPACE?



## NO? GOT CRAMPS?


Well no-body likes tight spaces so let us help you free some up in that utility center. (And score an awesome cookie with every order)

## SPACE SAVER SOLID BRASS MINI FILTER

- Save \$\$\$ by moving away from full size filters
- Solid Brass, Plated for esthetics
- Easy-access Filter disc
- 1/8" FPT or 10-32 Female ports
- Includes 35 Micron Filter disc, Sintered bronze
- Filter disc is readily available and inexpensive
- Disc measures 7/8" x 1/16"
- Better fit than full size filters in crowded junction boxes
- Mates nicely with our C-1836 Mini Regulator



V-1526	In-Line Space-Saver Filter, 1/8" FPT	\$18 MSRP
V-1530	In-Line Space-Saver Filter, 10-32 Female	\$18 MSRP
V-1527	35 Micron Replacement Disc	\$3 MSRP
V-1527-5	36 Micron Replacement Disc, Pkg of 5	\$14 MSRP



Delivering cookies is great, but remember . . . .

**We deliver Parts**

 **PARTS WAREHOUSE**  
*Let us be your dental resource*

1-800-817-6704  
[www.dentalparts.com](http://www.dentalparts.com)  
[CustomerService@Dentalparts.com](mailto:CustomerService@Dentalparts.com)